

Ziegelei 1 72336 Balingen-Frommern Germany

+0049-[0]7433-9933-0

- +0049-[0]7433-9933-149
- info@kern-sohn.com

Betriebsanleitung **KERN Schnittstellenadapter mit Kabel**für WLAN

KERN YKUP-05

Typ TYKUP-05-A Version 1.0 2021-09 D



Sie finden die aktuelle Version dieser Anleitung auch online unter: https://www.kern-sohn.com/shop/de/DOWNLOADS/ Unter der Rubrik Bedienungsanleitungen

TYKUP-05-A-BA-d-2110_WLAN



KERN Schnittstellenadapter mit Kabel Version 1.0 2021-09 Installationsanleitung für WLAN

Inhaltsverzeichnis

1	Allgemeines	2
1.1	Installation	2
1.2	Auszug aus KERN Communications Protocol KCP (Ref. manual 1.5.0)	5

1 Allgemeines

Kabellänge: 0,15 m

•1	 Über WLAN Es dürfen r schluss de 	können Wägedaten übertragen werden. Nur KERN KUP-Adapter an den 15-pol-Sub-D-An- r Waage angeschlossen werden!
	schluss de	r waage angeschlossen werden!

1.1 Installation

- Gerät ausschalten
- KUP-Adapter (WLAN) 15-pol-Sub-D-Anschluss des Gerätes einstecken
- Gerät einschalten
 - o Der KUP-Adapter wird automatisch vom Gerät erkannt
- ⇒ Nach dem Einschalten ohne Konfiguration erstellt das Gerät zunächst einen WLAN Accesspoint mit dem Namen "AI-Thinker_xxxxx.
- ⇒ Computer über diesen Accesspoint mit dem Gerät verbinden.
- ⇒ In einem Webbrowser die IP-Adresse 192.168.4.1 eingeben (Standard-IP). Es erscheint die Konfigurations-Website. Die statische IP wird über die KCP-Befehle vergeben.

Α	Betriebsmodus "apsta" auswählen
В	WLAN-Netzwerkname und entsprechendes Passwort eingeben
С	Einstellungen speichern und Ziel-Software neu starten (Reboot-Button)

(s. Abb.:)

ESP8266	WebConfig						Restore	Reb	oot
Serial Set	tting		SoftAP			Station			
Baud:	115200	~	SSID:	AI-THINKER_872B77		Mode:	apsta	~	Α
Databits:	8	~	Passwd:			AP Name:	PDWLAN		
Parity:	NONE	~	Auth Mode:	OPEN	~	AP Password:	12345678		Ь
Stopbits:	1	~	IP addr:	192.168.4.1		IP address:	0.0.0.0		
			Subnet mask:	255.255.255.0		Subnet mask:	0.0.0.0		
			Gateway:	192.168.4.1		Gateway:	0.0.0.0		
			Mac:	be:dd:c2:87:2b:77		Mac:	bc:dd:c2:87:2b:77		
		Save			Save		1	Save	С

D

Verbindung zum PC trennen und Stromversorgung des Geräts trennen.

1 Wurden Einstellungen am Gerät vorgenommen, ist darauf zu achten, dass das Gerät vollständig ausgeschalten ist! Die Einstellungen werden erst dann übernommen. Aktualisieren (Reboot-Button) und Speichern (Save-Button) sind nicht ausreichend!

Е	 Gerät wieder an die Stromversorgung anschließen,
	 Verbindung zum PC wieder herstellen,
	 Konfigurations-Webseite aufrufen und IP-Adresse pr üfen.

(s.Abb.:)

Serial Set	ting		SoftAP		Stati	on		
aud:	115200	~	SSID:	AI-THINKER_872B77	Mode	:	apsta	~
atabits:	8	~	Passwd:		AP Na	ame:	PDWLAN	
arity:	NONE	~	Auth Mode:	OPEN V	AP Pa	assword:	12345678	
topbits:	1	~	IP addr:	192.168.4.1	IP add	dress:	10.0.11.13	
			Subnet mask:	255.255.255.0	Subne	et mask:	255.255,0.0	
			Gateway:	192.168.4.1	Gatev	vay:	10.0.0.1	
			Mac:	be:dd:c2:87:2b:77	Mac:		bc:dd:c2:87:2b:77	
				_				

F	Konfigurations-Webseite schließen, PC mit dem ausgewählten Netzwerk verbinden
G	Ziel-Software (z.B. KERN Balance Connection) öffnen und IP-Adresse und Port 23 eintragen.

(s.Abb.:)

TCP/IP -> 10.0.11.13:23 < 440 (Gewichtswert-Parser) - Eigenschaften						
Allgemein Bus IP Port Eigenschaften TCP/UDP / IP Einstellungen:						
Verbindungsart: TCP - Client - verbindend Lokale IP Adresse: Port: Ziel Host/IP Adresse: 10.0.11.13 Keep-Alive: Image: Comparison of the second s						

Zur Abfrage der über DHCP zugewiesenen Konfiguration sowie zur spezifischen/statischen Konfiguration der IP-Adresse, Subnetz-Maske oder des Gateways können die KCP-Befehle **JNWx** verwendet werden.

1.2 Auszug aus KERN Communications Protocol KCP (Ref. manual 1.5.0)

JNWA - Query / set network address (IP) of WIFI Interface

Description

Use this command to query or set the network address (IP) of WIFI Interface.

Syntax

Command

JNWA	Query the current network address.
JNWA_«NetworkAddress»	Set the current network address.
JNWA_0.0.0	Activate DHCP.

Responses

JNWA_A_ <i>«NetworkAddress»</i>	Current network address (IP).
JNWAJA	Network address setting successfully performed.
JNWA_I	Command understood but currently not executable (device is currently executing another command, e.g. taring, or timeout as stability was not reached).
JNWA_L	Command understood but not executable (incor- rect parameter).

Parameters / Return values

Name	Туре	Values	Meaning
NetworkAddress	string		Network address (e.g. 192.168.0.1).

Comments

- All three commands, JNWA, JNWK and JNWG are required to enter sequentially for completing the setting of WIFI Interface.
- The exceptional case is activating the DHCP. The network mask and gateway address are not required. A single command "JNWA 0.0.0.0" can activate the DHCP of the WIFI Interface.
- It may take a few seconds to response to the command.

Examples

V	JNWA	Send current network address.
1	JNWA_A_192.168.0.1	The current network address is 192.168.0.1.
V	JNWA_192.168.0.1	Set network address to 192.168.0.1.
1	JNWA_A	Set network address setting successfully performed.
V	JNWA_0.0.0	Activate DHCP setting.
	JNWA_A	Successfully activated DHCP setting.

See also

→	JNWK - Query / set network mask
→	JNWG - Query / set gateway address

Description

Use this command to query or set the network mask of WIFI Interface.

Syntax

Command

JNWK	Query the current network mask.
JNWK _. «NetworkMask»	Set the current network mask.

Responses

JNWK_A_«NetworkMask»	Current network mask.
JNWK_A	Network mask setting successfully per-
	formed.
JNWKJI	Command understood but currently not
	executable
	(device is currently executing another
	command,
	e.g. taring, or timeout as stability was not
	reached).
JNWK_L	Command understood but not executable
	(incorrect
	parameter).

Parameters / Return values

Name	Туре	Values	Meaning
NetworkMask	string		Network mask (e.g. 255.255.255.0)

Comments

- All three commands, JNWA, JNWK and JNWG are required to enter sequentially for completing the setting of WIFI Interface.
- The exceptional case is activating the DHCP. The network mask and gateway address are not required. A single command "JNWA 0.0.0.0" can activate the DHCP of the WIFI Interface.
- It may take a few seconds to response to the command.

Examples

¥	JNWK	Send current network mask.
	JNWK_A_255.255.255.0	The current network mask is
		255.255.255.0.
\	JNWK_255.255.255.0	Set network mask to 255.255.255.0.
	JNWK_A	Set network mask setting successfully
		performed.

See also

→	JNWA - Query / set network address (IP)
>	JNWG - Query / set gateway address

Description

Use this command to query or set the gateway address of WIFI Interface.

Syntax

Command

JNWG	Query the current gateway address.
JNWG _. «GatewayAddress»	Set the current gateway address.

Responses

JNWG_A_«GatewayAddress»	Current gateway address.
JNWG_A	Gateway address setting successfully per-
	formed.
JNWG_I	Command understood but currently not
	executable
	(device is currently executing another
	command,
	e.g. taring, or timeout as stability was not
	reached).
JNWG_L	Command understood but not executable
	(incorrect
	parameter).

Parameters / Return values

Name	Туре	Values	Meaning
Gate-	string		Gateway address (e.g. 192.168.0.99)
wayAddress			

Comments

- All three commands, JNWA, JNWK and JNWG are required to enter sequentially for completing the setting of WIFI Interface.
- The exceptional case is activating the DHCP. The network mask and gateway address are not required. A single command "JNWA 0.0.0.0" can activate the DHCP of the WIFI Interface.
- It may take a few seconds to response to the command.

Examples

↓ ↓	JNWG	Send current gateway address.
1	JNWG_A_192.168.0.99	The current gateway address is 192.168.0.99.
Y	JNWG_192.168.0.99	Set gateway address to 192.168.0.99.
1	JNWG_A	Set gateway address setting successfully performed.

See also

→	JNWA - Query / set network address (IP)
>	JNWK - Query / set network mask